

Antibodies Product Datasheet

Product Name:	ACSL1 Antibody
Catalog Number:	C42475
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 210-290
Concentration:	1 mg/ml
Gene Name:	ACSL1 FACL1 FACL2 LACS LACS1 LACS2
Protein Name:	Long-chain-fatty-acid--CoA ligase 1 (EC 6.2.1.3) (Acyl-CoA synthetase 1) (ACS1) (Long-chain acyl-CoA synthetase 1) (LACS 1) (Long-chain acyl-CoA synthetase 2) (LACS 2) (Long-chain fatty acid-CoA ligase 2) (Palmitoyl-CoA ligase 1) (Palmitoyl-CoA ligase 2)
Other Name:	
Human Swiss Prot No.:	P33121
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	ACSL1 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year