

Antibodies Product Datasheet

Product Name:	ACTZ Antibody
Catalog Number:	C42484
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human ACTZ
Concentration:	1 mg/ml
Gene Name:	ACTR1A CTRN1
Protein Name:	Alpha-centractin (Centractin) (ARP1) (Actin-RPV) (Centrosome-associated actin homolog)
Other Name:	
Human Swiss Prot No.:	P61163
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human ACTZ
Storage and Stability:	-20°C/1 year