

Antibodies Product Datasheet

Product Name:	ADA11 Antibody
Catalog Number:	C42489
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human ADA11
Concentration:	1 mg/ml
Gene Name:	ADAM11 MDC
Protein Name:	ADA11
Other Name:	Disintegrin and metalloproteinase domain-containing protein 11 (ADAM 11) (Metalloproteinase-like, disintegrin-like, and cysteine-rich protein) (MDC)
Human Swiss Prot No.:	O75078
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.327% sodium azide.
Specificity:	This antibody detects endogenous levels of ADA11 at Human/Mouse
Storage and Stability:	-20°C/1 year