

Antibodies Product Datasheet

Product Name:	ADAM9 Antibody
Catalog Number:	C42501
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 120-200
Concentration:	1 mg/ml
Gene Name:	ADAM9 KIAA0021 MCMP MDC9 MLTNG
Protein Name:	Disintegrin and metalloproteinase domain-containing protein 9 (ADAM 9) (EC 3.4.24.-) (Cellular disintegrin-related protein) (Meltrin-gamma) (Metalloprotease/disintegrin/cysteine-rich protein 9) (Myeloma cell metalloproteinase)
Other Name:	
Human Swiss Prot No.:	Q13443
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	ADAM9 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year