

Antibodies Product Datasheet

Product Name:	ADCY2 Antibody
Catalog Number:	C42503
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 490-570
Concentration:	1 mg/ml
Gene Name:	ADCY2 KIAA1060
Protein Name:	Adenylate cyclase type 2 (EC 4.6.1.1) (ATP pyrophosphate-lyase 2) (Adenylate cyclase type II) (Adenylyl cyclase 2)
Other Name:	
Human Swiss Prot No.:	Q08462
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	ADCY2 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year