

## Antibodies Product Datasheet

|                               |                                                                                                                                                                                                                                                                       |
|-------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Name:</b>          | ADHX Antibody                                                                                                                                                                                                                                                         |
| <b>Catalog Number:</b>        | C42507                                                                                                                                                                                                                                                                |
| <b>Reactivity:</b>            | Human, Mouse, Rat                                                                                                                                                                                                                                                     |
| <b>Applications:</b>          | WB, ELISA                                                                                                                                                                                                                                                             |
| <b>Dilutions:</b>             | WB 1:500-2000 ELISA 1:5000-20000                                                                                                                                                                                                                                      |
| <b>Modification:</b>          | Unmodified/Total                                                                                                                                                                                                                                                      |
| <b>Source:</b>                | Polyclonal Rabbit                                                                                                                                                                                                                                                     |
| <b>Purification:</b>          | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.                                                                                                                                                 |
| <b>Immunogen:</b>             | Synthesized peptide derived from human protein . at AA range: 280-360                                                                                                                                                                                                 |
| <b>Concentration:</b>         | 1 mg/ml                                                                                                                                                                                                                                                               |
| <b>Gene Name:</b>             | ADH5 ADHX FDH                                                                                                                                                                                                                                                         |
| <b>Protein Name:</b>          | Alcohol dehydrogenase class-3 (EC 1.1.1.1) (Alcohol dehydrogenase 5) (Alcohol dehydrogenase class chi chain) (Alcohol dehydrogenase class-III) (Glutathione-dependent formaldehyde dehydrogenase) (FALDH) (FDH) (GSH-FDH) (EC 1.1.1.-) (S-(hydroxymethyl)glutathione) |
| <b>Other Name:</b>            |                                                                                                                                                                                                                                                                       |
| <b>Human Swiss Prot No.:</b>  | P11766                                                                                                                                                                                                                                                                |
| <b>Formulation:</b>           | Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.                                                                                                                                                                                                        |
| <b>Specificity:</b>           | ADHX Polyclonal Antibody detects endogenous levels of protein.                                                                                                                                                                                                        |
| <b>Storage and Stability:</b> | -20°C/1 year                                                                                                                                                                                                                                                          |