

Antibodies Product Datasheet

Product Name:	ADT2 Antibody
Catalog Number:	C42509
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 40-120
Concentration:	1 mg/ml
Gene Name:	SLC25A5 ANT2
Protein Name:	ADP/ATP translocase 2 (ADP,ATP carrier protein 2) (ADP,ATP carrier protein, fibroblast isoform) (Adenine nucleotide translocator 2) (ANT 2) (Solute carrier family 25 member 5)
Other Name:	
Human Swiss Prot No.:	P05141
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	ADT2 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year