

Antibodies Product Datasheet

Product Name:	AF1L1 Antibody
Catalog Number:	C42513
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human AF1L1
Concentration:	1 mg/ml
Gene Name:	AFAP1L1
Protein Name:	AF1L1
Other Name:	Actin filament-associated protein 1-like 1 (AFAP1-like protein 1)
Human Swiss Prot No.:	Q8TED9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.325% sodium azide.
Specificity:	This antibody detects endogenous levels of AF1L1 at Human/Mouse/Rat
Storage and Stability:	-20°C/1 year