

Antibodies Product Datasheet

Product Name:	AFF2 Antibody
Catalog Number:	C42517
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human AFF2
Concentration:	1 mg/ml
Gene Name:	AFF2 FMR2 OX19
Protein Name:	AF4/FMR2 family member 2 (Fragile X E mental retardation syndrome protein) (Fragile X mental retardation 2 protein) (FMR2P) (Protein FMR-2) (Protein Ox19)
Other Name:	
Human Swiss Prot No.:	P51816
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human AFF2
Storage and Stability:	-20°C/1 year