

Antibodies Product Datasheet

Product Name:	AGFG1 Antibody
Catalog Number:	C42521
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	AGFG1 HRB RAB RIP
Protein Name:	Arf-GAP domain and FG repeat-containing protein 1 (HIV-1 Rev-binding protein) (Nucleoporin-like protein RIP) (Rev-interacting protein) (Rev/Rex activation domain-binding protein)
Other Name:	
Human Swiss Prot No.:	P52594
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	AGFG1 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year