

Antibodies Product Datasheet

Product Name:	AGK Antibody
Catalog Number:	C42522
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 330-410
Concentration:	1 mg/ml
Gene Name:	AGK MULK
Protein Name:	Acylglycerol kinase, mitochondrial (hAGK) (EC 2.7.1.107) (EC 2.7.1.94) (Multiple substrate lipid kinase) (HsMuLK) (MuLK) (Multi-substrate lipid kinase)
Other Name:	
Human Swiss Prot No.:	Q53H12
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	AGK Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year