

Antibodies Product Datasheet

Product Name:	AGT2 Antibody
Catalog Number:	C42526
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 340-420
Concentration:	1 mg/ml
Gene Name:	AGXT2 AGT2
Protein Name:	Alanine--glyoxylate aminotransferase 2, mitochondrial (AGT 2) (EC 2.6.1.44) ((R)-3-amino-2-methylpropionate--pyruvate transaminase) (EC 2.6.1.40) (Beta-ALAAT II) (Beta-alanine-pyruvate aminotransferase) (D-AIBAT)
Other Name:	
Human Swiss Prot No.:	Q9BYV1
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	AGT2 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year