

Antibodies Product Datasheet

Product Name:	AK1C4 Antibody
Catalog Number:	C42532
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 1-80
Concentration:	1 mg/ml
Gene Name:	AKR1C4 CHDR
Protein Name:	Aldo-keto reductase family 1 member C4 (EC 1.1.1.-) (3-alpha-HSD1) (3-alpha-hydroxysteroid dehydrogenase type I) (EC 1.1.1.50) (Chlordecone reductase) (CDR) (EC 1.1.1.225) (Dihydrodiol dehydrogenase 4) (DD-4) (DD4) (HAKRA)
Other Name:	
Human Swiss Prot No.:	P17516
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	AK1C4 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year