

Antibodies Product Datasheet

Product Name:	AKAP4 Antibody
Catalog Number:	C42533
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 430-510
Concentration:	1 mg/ml
Gene Name:	AKAP4 AKAP82
Protein Name:	A-kinase anchor protein 4 (AKAP-4) (A-kinase anchor protein 82 kDa) (AKAP 82) (hAKAP82) (Major sperm fibrous sheath protein) (HI) (Protein kinase A-anchoring protein 4) (PRKA4)
Other Name:	
Human Swiss Prot No.:	Q5JQC9
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	AKAP4 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year