

Antibodies Product Datasheet

Product Name:	AKAP6 Antibody
Catalog Number:	C42534
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section)
Dilutions:	IHC-p 1:50-300
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 30-110
Concentration:	1 mg/ml
Gene Name:	AKAP6 AKAP100 KIAA0311
Protein Name:	A-kinase anchor protein 6 (AKAP-6) (A-kinase anchor protein 100 kDa) (AKAP 100) (Protein kinase A-anchoring protein 6) (PRKA6) (mAKAP)
Other Name:	
Human Swiss Prot No.:	Q13023
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	AKAP6 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year