

Antibodies Product Datasheet

Product Name:	AKAP9 Antibody
Catalog Number:	C42535
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section)
Dilutions:	IHC-p 1:50-300
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 2610-2690
Concentration:	1 mg/ml
Gene Name:	AKAP9 AKAP350 AKAP450 KIAA0803
Protein Name:	A-kinase anchor protein 9 (AKAP-9) (A-kinase anchor protein 350 kDa) (AKAP 350) (hgAKAP 350) (A-kinase anchor protein 450 kDa) (AKAP 450) (AKAP 120-like protein) (Centrosome- and Golgi-localized PKN-associated protein) (CG-NAP) (Protein hyperion) (Protein
Other Name:	
Human Swiss Prot No.:	Q99996
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	AKAP9 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year