

Antibodies Product Datasheet

Product Name:	AL2CL Antibody
Catalog Number:	C42540
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human AL2CL
Concentration:	1 mg/ml
Gene Name:	ALS2CL
Protein Name:	AL2CL
Other Name:	ALS2 C-terminal-like protein
Human Swiss Prot No.:	Q60I27
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.222% sodium azide.
Specificity:	This antibody detects endogenous levels of AL2CL at Human/Mouse
Storage and Stability:	-20°C/1 year