

Antibodies Product Datasheet

Product Name:	ALKB8 Antibody
Catalog Number:	C42547
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human ALKB8
Concentration:	1 mg/ml
Gene Name:	ALKBH8 ABH8
Protein Name:	Alkylated DNA repair protein alkB homolog 8 (EC 1.14.11.-) (Probable alpha-ketoglutarate-dependent dioxygenase ABH8) (S-adenosyl-L-methionine-dependent tRNA methyltransferase ABH8) (tRNA (carboxymethyluridine(34)-5-O)-methyltransferase ABH8) (EC 2.1.1.229
Other Name:	
Human Swiss Prot No.:	Q96BT7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human ALKB8
Storage and Stability:	-20°C/1 year