

Antibodies Product Datasheet

Product Name:	ANO4 Antibody
Catalog Number:	C42577
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human ANO4
Concentration:	1 mg/ml
Gene Name:	ANO4 TMEM16D
Protein Name:	ANO4
Other Name:	Anoctamin-4 (Transmembrane protein 16D)
Human Swiss Prot No.:	Q32M45
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.232% sodium azide.
Specificity:	This antibody detects endogenous levels of ANO4 at Human/Mouse
Storage and Stability:	-20°C/1 year