

Antibodies Product Datasheet

Product Name:	ANR28 Antibody
Catalog Number:	C42582
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 980-1060
Concentration:	1 mg/ml
Gene Name:	ANKRD28 KIAA0379
Protein Name:	Serine/threonine-protein phosphatase 6 regulatory ankyrin repeat subunit A (PP6-ARS-A) (Serine/threonine-protein phosphatase 6 regulatory subunit ARS-A) (Ankyrin repeat domain-containing protein 28) (Phosphatase interactor targeting protein hnRNP K) (PITK)
Other Name:	
Human Swiss Prot No.:	O15084
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	ANR28 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year