

Antibodies Product Datasheet

Product Name:	ANXA4 Antibody
Catalog Number:	C42593
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	ANXA4 ANX4
Protein Name:	Annexin A4 (35-beta calcimedlin) (Annexin IV) (Annexin-4) (Carbohydrate-binding protein p33/p41) (Chromobindin-4) (Endonexin I) (Lipocortin IV) (P32.5) (PP4-X) (Placental anticoagulant protein II) (PAP-II) (Protein II)
Other Name:	
Human Swiss Prot No.:	P09525
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	ANXA4 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year