

Antibodies Product Datasheet

Product Name:	AOFB Antibody
Catalog Number:	C42596
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: Internal
Concentration:	1 mg/ml
Gene Name:	MAOB
Protein Name:	Amine oxidase [flavin-containing] B (EC 1.4.3.4) (Monoamine oxidase type B) (MAO-B)
Other Name:	
Human Swiss Prot No.:	P27338
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	AOFB Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year