

Antibodies Product Datasheet

Product Name:	AP1B1 Antibody
Catalog Number:	C42597
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human AP1B1
Concentration:	1 mg/ml
Gene Name:	AP1B1 ADTB1 BAM22 CLAPB2
Protein Name:	AP1B1
Other Name:	AP-1 complex subunit beta-1 (Adapter-related protein complex 1 subunit beta-1) (Adaptor protein complex AP-1 subunit beta-1) (Beta-1-adaptin) (Beta-adaptin 1) (Clathrin assembly protein complex 1 beta large chain) (Golgi adaptor HA1/AP1 adaptin beta subun
Human Swiss Prot No.:	Q10567
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.201% sodium azide.
Specificity:	This antibody detects endogenous levels of AP1B1 at Human/Mouse/Rat
Storage and Stability:	-20°C/1 year