

Antibodies Product Datasheet

Product Name:	AP1G1 Antibody
Catalog Number:	C42598
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	AP1G1 ADTG CLAPG1
Protein Name:	AP-1 complex subunit gamma-1 (Adapter-related protein complex 1 subunit gamma-1) (Adaptor protein complex AP-1 subunit gamma-1) (Clathrin assembly protein complex 1 gamma-1 large chain) (Gamma1-adaptin) (Golgi adaptor HA1/AP1 adaptin subunit gamma-1)
Other Name:	
Human Swiss Prot No.:	O43747
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	AP1G1 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year