

Antibodies Product Datasheet

Product Name:	AP2D Antibody
Catalog Number:	C42602
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 280-360
Concentration:	1 mg/ml
Gene Name:	TFAP2D TFAP2BL1
Protein Name:	Transcription factor AP-2-delta (AP2-delta) (Activating enhancer-binding protein 2-delta) (Transcription factor AP-2-beta-like 1)
Other Name:	
Human Swiss Prot No.:	Q7Z6R9
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	AP2D Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year