

Antibodies Product Datasheet

Product Name:	AP2M1 Antibody
Catalog Number:	C42603
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 100-180
Concentration:	1 mg/ml
Gene Name:	AP2M1 CLAPM1 KIAA0109
Protein Name:	AP-2 complex subunit mu (AP-2 mu chain) (Adapter-related protein complex 2 mu subunit) (Adaptin-mu2) (Adaptor protein complex AP-2 subunit mu) (Clathrin assembly protein complex 2 medium chain) (Clathrin coat assembly protein AP50) (Clathrin coat-associat
Other Name:	
Human Swiss Prot No.:	Q96CW1
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	AP2M1 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year