

Antibodies Product Datasheet

Product Name:	APBA1 Antibody
Catalog Number:	C42604
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 180-260
Concentration:	1 mg/ml
Gene Name:	APBA1 MINT1 X11
Protein Name:	Amyloid beta A4 precursor protein-binding family A member 1 (Adapter protein X11alpha) (Neuron-specific X11 protein) (Neuronal Munc18-1-interacting protein 1) (Mint-1)
Other Name:	
Human Swiss Prot No.:	Q02410
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	APBA1 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year