

Antibodies Product Datasheet

Product Name:	APOH Antibody
Catalog Number:	C42614
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	APOH B2G1
Protein Name:	Beta-2-glycoprotein 1 (APC inhibitor) (Activated protein C-binding protein) (Anticardiolipin cofactor) (Apolipoprotein H) (Apo-H) (Beta-2-glycoprotein I) (B2GPI) (Beta(2)GPI)
Other Name:	
Human Swiss Prot No.:	P02749
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	APOH Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year