

Antibodies Product Datasheet

Product Name:	ARAP1 Antibody
Catalog Number:	C42621
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	ARAP1 CENTD2 KIAA0782
Protein Name:	Arf-GAP with Rho-GAP domain, ANK repeat and PH domain-containing protein 1 (Centaurin-delta-2) (Cnt-d2)
Other Name:	
Human Swiss Prot No.:	Q96P48
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	ARAP1 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year