

Antibodies Product Datasheet

Product Name:	ARHG6 Antibody
Catalog Number:	C42628
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 580-660
Concentration:	1 mg/ml
Gene Name:	ARHGEF6 COOL2 KIAA0006 PIXA
Protein Name:	Rho guanine nucleotide exchange factor 6 (Alpha-Pix) (COOL-2) (PAK-interacting exchange factor alpha) (Rac/Cdc42 guanine nucleotide exchange factor 6)
Other Name:	
Human Swiss Prot No.:	Q15052
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	ARHG6 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year