

Antibodies Product Datasheet

Product Name:	ARI1A Antibody
Catalog Number:	C42632
Reactivity:	Human, Mouse
Applications:	IHC-p, IF(paraffin section)
Dilutions:	IHC-p 1:50-300
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 1270-1350
Concentration:	1 mg/ml
Gene Name:	ARID1A BAF250 BAF250A C1orf4 OSA1 SMARCF1
Protein Name:	AT-rich interactive domain-containing protein 1A (ARID domain-containing protein 1A) (B120) (BRG1-associated factor 250) (BAF250) (BRG1-associated factor 250a) (BAF250A) (Osa homolog 1) (hOSA1) (SWI-like protein) (SWI/SNF complex protein p270) (SWI/SNF-re
Other Name:	
Human Swiss Prot No.:	O14497
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	ARI1A Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year