

Antibodies Product Datasheet

Product Name:	ARI3A Antibody
Catalog Number:	C42633
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	ARID3A DRIL1 DRIL3 DRX E2FBP1
Protein Name:	AT-rich interactive domain-containing protein 3A (ARID domain-containing protein 3A) (B-cell regulator of IgH transcription) (Bright) (Dead ringer-like protein 1) (E2F-binding protein 1)
Other Name:	
Human Swiss Prot No.:	Q99856
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	ARI3A Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year