

Antibodies Product Datasheet

Product Name:	ARI3B Antibody
Catalog Number:	C42634
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human ARI3B
Concentration:	1 mg/ml
Gene Name:	ARID3B BDP DRIL2
Protein Name:	ARI3B
Other Name:	AT-rich interactive domain-containing protein 3B (ARID domain-containing protein 3B) (Bright and dead ringer protein) (Bright-like protein)
Human Swiss Prot No.:	Q8IVW6
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.83% sodium azide.
Specificity:	This antibody detects endogenous levels of ARI3B at Human/Mouse
Storage and Stability:	-20°C/1 year