

Antibodies Product Datasheet

Product Name:	ARL14 Antibody
Catalog Number:	C42636
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human ARL14
Concentration:	1 mg/ml
Gene Name:	ARL14 ARF7
Protein Name:	ADP-ribosylation factor-like protein 14 (ADP-ribosylation factor 7)
Other Name:	
Human Swiss Prot No.:	Q8N4G2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human ARL14
Storage and Stability:	-20°C/1 year