

Antibodies Product Datasheet

Product Name:	ARPC2 Antibody
Catalog Number:	C42643
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human ARPC2
Concentration:	1 mg/ml
Gene Name:	ARPC2 ARC34 PRO2446
Protein Name:	Actin-related protein 2/3 complex subunit 2 (Arp2/3 complex 34 kDa subunit) (p34-ARC)
Other Name:	
Human Swiss Prot No.:	O15144
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human ARPC2
Storage and Stability:	-20°C/1 year