

Antibodies Product Datasheet

Product Name:	ASAP2 Antibody
Catalog Number:	C42650
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 640-720
Concentration:	1 mg/ml
Gene Name:	ASAP2 DDEF2 KIAA0400
Protein Name:	Arf-GAP with SH3 domain, ANK repeat and PH domain-containing protein 2 (Development and differentiation-enhancing factor 2) (Paxillin-associated protein with ARF GAP activity 3) (PAG3) (Pyk2 C-terminus-associated protein) (PAP)
Other Name:	
Human Swiss Prot No.:	O43150
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	ASAP2 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year