

Antibodies Product Datasheet

Product Name:	ASCL4 Antibody
Catalog Number:	C42658
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human ASCL4
Concentration:	1 mg/ml
Gene Name:	ASCL4 BHLHA44 HASH4
Protein Name:	Achaete-scute homolog 4 (ASH-4) (hASH4) (Achaete-scute-like protein 4) (Class A basic helix-loop-helix protein 44) (bHLHa44)
Other Name:	
Human Swiss Prot No.:	Q6XD76
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human ASCL4
Storage and Stability:	-20°C/1 year