

Antibodies Product Datasheet

Product Name:	ASH1L Antibody
Catalog Number:	C42660
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human ASH1L
Concentration:	1 mg/ml
Gene Name:	ASH1L KIAA1420 KMT2H
Protein Name:	Histone-lysine N-methyltransferase ASH1L (EC 2.1.1.43) (ASH1-like protein) (huASH1) (Absent small and homeotic disks protein 1 homolog) (Lysine N-methyltransferase 2H)
Other Name:	
Human Swiss Prot No.:	Q9NR48
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human ASH1L
Storage and Stability:	-20°C/1 year