

Antibodies Product Datasheet

Product Name:	ASNS Antibody
Catalog Number:	C42664
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human ASNS
Concentration:	1 mg/ml
Gene Name:	ASNS TS11
Protein Name:	ASNS
Other Name:	Asparagine synthetase [glutamine-hydrolyzing] (EC 6.3.5.4) (Cell cycle control protein TS11) (Glutamine-dependent asparagine synthetase)
Human Swiss Prot No.:	P08243
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.84% sodium azide.
Specificity:	This antibody detects endogenous levels of ASNS at Human/Mouse/Rat
Storage and Stability:	-20°C/1 year