

Antibodies Product Datasheet

Product Name:	ASPP2 Antibody
Catalog Number:	C42666
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 60-140
Concentration:	1 mg/ml
Gene Name:	TP53BP2 ASPP2 BBP
Protein Name:	Apoptosis-stimulating of p53 protein 2 (Bcl2-binding protein) (Bbp) (Renal carcinoma antigen NY-REN-51) (Tumor suppressor p53-binding protein 2) (53BP2) (p53-binding protein 2) (p53BP2)
Other Name:	
Human Swiss Prot No.:	Q13625
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	ASPP2 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year