

Antibodies Product Datasheet

Product Name:	AT2B2 Antibody
Catalog Number:	C42673
Reactivity:	Human, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 250-330
Concentration:	1 mg/ml
Gene Name:	ATP2B2 PMCA2
Protein Name:	Plasma membrane calcium-transporting ATPase 2 (PMCA2) (EC 3.6.3.8) (Plasma membrane calcium ATPase isoform 2) (Plasma membrane calcium pump isoform 2)
Other Name:	
Human Swiss Prot No.:	Q01814
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	AT2B2 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year