

Antibodies Product Datasheet

Product Name:	AT2B3 Antibody
Catalog Number:	C42674
Reactivity:	H, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 140-220
Concentration:	1 mg/ml
Gene Name:	ATP2B3
Protein Name:	Plasma membrane calcium-transporting ATPase 3 (PMCA3) (EC 3.6.3.8) (Plasma membrane calcium ATPase isoform 3) (Plasma membrane calcium pump isoform 3)
Other Name:	
Human Swiss Prot No.:	Q16720
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	AT2B3 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year