

Antibodies Product Datasheet

Product Name:	AT8B2 Antibody
Catalog Number:	C42681
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human AT8B2
Concentration:	1 mg/ml
Gene Name:	ATP8B2 ATP10A KIAA1137
Protein Name:	Probable phospholipid-transporting ATPase ID (EC 3.6.3.1) (ATPase class I type 8B member 2)
Other Name:	
Human Swiss Prot No.:	P98198
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human AT8B2
Storage and Stability:	-20°C/1 year