

Antibodies Product Datasheet

Product Name:	ATL2 Antibody
Catalog Number:	C42685
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human ATL2
Concentration:	1 mg/ml
Gene Name:	ADAMTSL2 KIAA0605
Protein Name:	ATL2
Other Name:	ADAMTS-like protein 2 (ADAMTSL-2)
Human Swiss Prot No.:	Q86TH1
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.216% sodium azide.
Specificity:	This antibody detects endogenous levels of ATL2 at Human/Mouse
Storage and Stability:	-20°C/1 year