

Antibodies Product Datasheet

Product Name:	ATOH1 Antibody
Catalog Number:	C42687
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human ATOH1
Concentration:	1 mg/ml
Gene Name:	ATOH1 ATH1 BHLHA14
Protein Name:	ATOH1
Other Name:	Protein atonal homolog 1 (Class A basic helix-loop-helix protein 14) (bHLHa14) (Helix-loop-helix protein hATH-1) (hATH1)
Human Swiss Prot No.:	Q92858
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.46% sodium azide.
Specificity:	This antibody detects endogenous levels of ATOH1 at Human/Mouse
Storage and Stability:	-20°C/1 year