

## Antibodies Product Datasheet

|                               |                                                                                                                          |
|-------------------------------|--------------------------------------------------------------------------------------------------------------------------|
| <b>Product Name:</b>          | ATOH7 Antibody                                                                                                           |
| <b>Catalog Number:</b>        | C42688                                                                                                                   |
| <b>Reactivity:</b>            | Human, Mouse                                                                                                             |
| <b>Applications:</b>          | WB, ELISA                                                                                                                |
| <b>Dilutions:</b>             | WB 1:500-2000 ELISA 1:5000-20000                                                                                         |
| <b>Modification:</b>          | Unmodified/Total                                                                                                         |
| <b>Source:</b>                | Polyclonal Rabbit                                                                                                        |
| <b>Purification:</b>          | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.    |
| <b>Immunogen:</b>             | Synthesized peptide derived from part region of human protein                                                            |
| <b>Concentration:</b>         | 1 mg/ml                                                                                                                  |
| <b>Gene Name:</b>             | ATOH7 ATH5 BHLHA13                                                                                                       |
| <b>Protein Name:</b>          | Protein atonal homolog 7 (Class A basic helix-loop-helix protein 13) (bHLHa13) (Helix-loop-helix protein hATH-5) (hATH5) |
| <b>Other Name:</b>            |                                                                                                                          |
| <b>Human Swiss Prot No.:</b>  | Q8N100                                                                                                                   |
| <b>Formulation:</b>           | Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.                                                           |
| <b>Specificity:</b>           | ATOH7 Polyclonal Antibody detects endogenous levels of protein.                                                          |
| <b>Storage and Stability:</b> | -20°C/1 year                                                                                                             |