

Antibodies Product Datasheet

Product Name:	ATS13 Antibody
Catalog Number:	C42697
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 940-1020
Concentration:	1 mg/ml
Gene Name:	ADAMTS13 C9orf8 UNQ6102/PRO20085
Protein Name:	A disintegrin and metalloproteinase with thrombospondin motifs 13 (ADAM-TS 13) (ADAM-TS13) (ADAMTS-13) (EC 3.4.24.87) (von Willebrand factor-cleaving protease) (vWF-CP) (vWF-cleaving protease)
Other Name:	
Human Swiss Prot No.:	Q76LX8
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	ATS13 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year