

Antibodies Product Datasheet

Product Name:	ATS15 Antibody
Catalog Number:	C42699
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human ATS15
Concentration:	1 mg/ml
Gene Name:	ADAMTS15
Protein Name:	ATS15
Other Name:	A disintegrin and metalloproteinase with thrombospondin motifs 15 (ADAM-TS 15) (ADAM-TS15) (ADAMTS-15) (EC 3.4.24.-)
Human Swiss Prot No.:	Q8TE58
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.207% sodium azide.
Specificity:	This antibody detects endogenous levels of ATS15 at Human/Mouse
Storage and Stability:	-20°C/1 year