

Antibodies Product Datasheet

Product Name:	ATS5 Antibody
Catalog Number:	C42703
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	ADAMTS5 ADAMTS11 ADMP2
Protein Name:	A disintegrin and metalloproteinase with thrombospondin motifs 5 (ADAM-TS 5) (ADAM-TS5) (ADAMTS-5) (EC 3.4.24.-) (A disintegrin and metalloproteinase with thrombospondin motifs 11) (ADAM-TS 11) (ADAMTS-11) (ADMP-2) (Aggrecanase-2)
Other Name:	
Human Swiss Prot No.:	Q9UNA0
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	ATS5 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year