

Antibodies Product Datasheet

Product Name:	ATX3L Antibody
Catalog Number:	C42707
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human ATX3L
Concentration:	1 mg/ml
Gene Name:	ATXN3L ATX3L MJDL
Protein Name:	ATX3L
Other Name:	Putative ataxin-3-like protein (EC 3.4.19.12) (Machado-Joseph disease protein 1-like)
Human Swiss Prot No.:	Q9H3M9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.64% sodium azide.
Specificity:	This antibody detects endogenous levels of ATX3L at Human
Storage and Stability:	-20°C/1 year